

Carestream Health, Inc.

Exclusive Manufacturer of **Kodak** Dental Systems

Media Contact:

Emma Gyllström
+49 (0)711 20707-586
emma.gyllstroem@cshdental.com

New KODAK 9500 Cone Beam 3D System Redefines CBCT Design, Ease-of-Use

Medium, Large Field 3D Programs Deliver Dual Jaw, Craniofacial Imaging

COLOGNE, March 24 — The KODAK 9500 Cone Beam 3D System was introduced by Carestream Health, Inc., for easy-to-use three-dimensional imaging of the dental and facial anatomy. Built on the proven and reliable KODAK 9000 platform, the KODAK 9500 3D System enables dental professionals to quickly produce magnificent 3D images, ranging from dual jaw to full craniofacial images, while optimizing dose.

The KODAK 9500 3D System is well suited for general dentists and dental specialists performing orthodontic, periodontic and oral and maxillofacial surgical (OMS) procedures including implants, tooth extractions and surgical planning. With high-quality, anatomically-correct 3D images down to 0.2 mm voxel (VOLUME piXEL) size, dental professionals are able to leverage the incredible image detail of 3D to provide more accurate diagnoses, improved treatment planning and more efficient patient care.

“The KODAK 9500 3D system is among the most user- and patient-friendly large-field CBCT systems I’ve seen, said Dr. Dirk Schulze, Director of the Department of Oral Radiology, Department of Oral and Maxillofacial Surgery, University Hospital Freiburg (Germany). “My team really appreciates its user-friendliness and the large size of the volume.”

Smart Design Enhances Image Quality, Patient Experience

Featuring a user-friendly and stylish design that aids in patient comfort and satisfaction, the KODAK 9500 3D System provides innovative features like face-to-face patient positioning to accommodate standing or seated patients, including wheelchair patients, of varying body types and sizes. Program selection with pre-programmed

settings occurs directly on the system's intuitive interface, simplifying patient setup and exam preparation.

The system also provides full-featured 3D visualization software designed for the needs of dental professionals. The software is intuitive and easy-to-use, enabling users to focus on diagnosis and treatment. The application includes an implant planning/simulation module and implant library, and will be compatible with third-party implant planning applications. Image output occurs in DICOM formats, supporting flexible output options including printing and customizable report templates.

By maximizing ease-of-use and patient comfort, users can get a high-quality 3D image with the first take.

Anatomical Programs Optimize Patient Care, Safety

The KODAK 9500 3D System captures dual jaw (9 cm x 15 cm) or all dento-maxillofacial anatomy (18.4 cm x 20.6 cm) in one acquisition for a wide range of clinical applications. With the additional ability to control patient dose through variable settings of mA and kV, dental professionals can manage the patient dose by focusing on their immediate clinical needs. This enables users to maximize image detail required for the specific application while managing patient dose.

Pricing and Availability

The KODAK 9500 3D System is currently available for purchase in limited quantities in parts of Europe, with general availability scheduled in May of this year. For more information on pricing and availability, as well as a complete description of system features and benefits, contact your Sales and Service Representative at 00800 4567 7654 or go to www.kodakdental.com/9500.

###

About KODAK Dental Systems and Carestream Health, Inc.

Carestream Health, Inc., a leading provider of dental and medical imaging and IT solutions, is the exclusive maker of KODAK Dental Systems. KODAK Dental Imaging and Practice Management Systems help dental professionals streamline workflow and enhance patient care by enabling them to capture, share, and use images and information easily and effectively.

Headquartered in Rochester, N.Y., and serving customers in more than 150 countries, Carestream Health is an independent company operating within the family of companies owned by Onex Corporation (TSX: OCX.TO) of Toronto, Canada.

For more information, contact your Carestream Health representative or visit www.kodakdental.com.